

## Diagnostic Test

This invention relates to a cellular diagnostic test for the detection of infection and inflammation.

5

It is difficult to name a disease where inflammation is not present. The main actor in the creation and perpetuation of inflammation is the blood white cell (leucocyte) named polymorphonuclear neutrophil (PMN). PMNs are by far the most numerous of blood leucocytes. They also play a key role in the defence against microbes, as shown by the onset of deadly infections when they are impaired or too low in number. The paradox is that for this pivotal cell there exists no biological test that is both recognised by academic medicine and routinely accessible to the patients.

Consequently there are many gaps when it comes to understand the pathogenesis of many diseases, in turn resulting in the lack of proper diagnostic procedures. A good example can be found in a syndrome named "food intolerance" or "food sensitivity" that describes any clinically abnormal response attributed to an exposure to a food or food component, whether immunological or non-immunological.

Food intolerance reflecting PMN activation may be an underlying mechanism contributing to some serious chronic and even autoimmune diseases. Furthermore, several studies have shown that many conditions such as migraine, irritable bowel syndrome, rheumatoid arthritis, respond to dietary modification. However, the offending foods are not identifiable by usual immunological methods.

Several methods of assaying leucocyte activation in whole blood are offered to selected groups of patents, generally restricted to food intolerance:

1) the antigen leucocyte cellular antibody test (ALCAT), which is intended for the determination of the effects of multipathogenic mechanisms involved in food sensitivities. This method suffers from the disadvantage that it requires costly and specialised equipment and relatively large samples of whole blood, e.g. between 15 and 20 mL.

2) The Yorke Test is confined to measuring a class of antibodies, the IgG, based on the debatable assumption that food intolerance is exclusively a reaction due to the immune system.

3) Other methods measure the release of chemical mediators. One such example is the Food Allergen Cellular Test (FACT), where food sensitised leucocytes release an inflammatory group of molecules, the leukotrienes, which it is claimed reflect allergic responses but are unable to distinguish between food allergy and intolerance.

Practically all of the prior art methods of diagnosing PMN activation, e.g. in food intolerance are indirect methods which do not address directly the activation of the cell.

By contrast, we have invented and developed a new diagnostic method, which overcomes or substantially mitigates the disadvantages of prior art methods, whilst at the same time giving reliable measures of the state of PMN. It then becomes possible to explore infection,

inflammation and other syndromes such as food intolerance by directly measuring PMN activation.

According to the invention we provide a method of  
5 diagnosing the level of PMN activation following incubation of a leucocyte sample prepared from whole blood. PMN activation may be either spontaneous (thus reflecting the basic state of the cell in circulating blood) or after exposure to a classical PMN activating  
10 agent such as phorbol-myristate acetate (PMA) or to substances suspected of acting on PMNs such as food extracts.

The method allows detecting whether or not PMNs in the sample have become activated following incubation in  
15 the presence of the said compounds.

The invention relates to a method of diagnosing intolerance by a subject to a specified substance which comprises incubating a leucocyte sample prepared from  
20 whole blood drawn from the subject with an extract of said substance and detecting whether or not PMNs in the sample have become activated.

A number of methods can be used to determine if PMNs have become activated.

25 The change in PMN morphology which results from activation can be determined optically, e.g. using a microscope.

According to an embodiment, the PMN activation is detected optically.

30 However, although this method is accurate, it is very time consuming, which makes it unfit for routine testing. Alternative methods of detecting PMN activation include luminescence, detecting the uptake of dyes or fluorescent

markers, photometry, detecting the release of cytokines or other bioactive molecules (e.g. leukotrienes), microbial proteins (e.g. defensins) and/or free radicals (e.g. reactive oxygen metabolites).

5       According to an embodiment, PMN activation is detected by determining whether they adhere to surface. This makes use of the fact that PMNs tend to become more adhesive on activation. As such, a preferred method of detecting PMN activation is to monitor the increased  
10       adhesiveness of activated PMNs, relative to inactivated ones. In particular, activated PMNs tend to adhere to plastic surfaces, e.g. the surface of a well in a 96-well titre plate or in a sample tube, for example an Eppendorf® tube. Thus, PMN adhesion may be determined by assaying  
15       adhesion to a plastic surface that may be a plastic multi-well titre plate that may contains 96 wells. Accordingly, the number of PMNs adhering to a plastic surface can be used as a measure of PMN activation.

Adhesion may also be detected by lysis and assaying  
20       for one or more intracellular markers for PMNs. The one or more markers may be selected from acid hydrolases, myeloperoxidases, lysozyme, lactoferrin, neutral proteases, serine proteases and lactic dehydrogenase.

According to an embodiment, the plastic surface is  
25       washed to remove unreacted cells, after incubation of the sample with the food extract but before assaying adhesion.

In fact, in general, when the PMN adhesiveness is determined by detecting adhesion to any plastic surface, such as the well of 96-well plate, non-adhering cells are  
30       removed by vigorous washing. The number of PMNs adhering to the plastic surface can be assayed by several methods, e.g. those methods mentioned above.

When the surface is a well, e.g. in a 96-well titre plate, the well can be scored turbidometrically. However, preferably, the number of cells remaining after adhesion and washing is detected by lysis and assaying for one or more intracellular PMN components.

Methods of cell lysis are well known in the art and include adding a detergent, e.g. Triton X100, to adhered cells, thereby releasing intracellular components into the cell suspending medium. Several intracellular PMN markers for activation include acid hydrolases, myeloperoxidases, lysozyme, lactoferrin, neutral and serine proteases. A preferred marker is the cytoplasmic enzyme lactic dehydrogenase (LDH), as the ratio of extracellular to intracellular enzyme is high, of the order 1:500.

15

The method is preferably carried out using a 96 well titre plate, allowing an array of food extracts to be investigated in run, together with suitable controls.

## 20 **Example**

The following example is about detection of neutrophil reactivity to foodstuffs. All substances capable of activating neutrophils e.g. PMA, (always used as a positive control, see below), anaphylatoxins, lectins, hetero- or auto-antigens... can be used in the same manner as food stuffs to measure their effect on neutrophil adherence.

### 30 **1. Preparing semi-purified leucocytes**

A sample of whole blood (2 to 10 ml) is drawn from subjects being tested by phlebotomy into plastic tubes

containing a  $\text{Ca}^{2+}$ -chelating anticoagulant e.g. CPDA or ACD. An aliquot of 270 $\mu\text{L}$  Plasmagel® is added followed by gently mixing. Tubes are tilted at an angle of  $45^\circ$  and left for 1 hour in this position. The erythrocytes (red blood cells) sediment during this time leaving the white blood cell fraction in the plasma/plasmagel layer which now constitutes the supernatant.

Cells in supernatant are counted using a microscope: a 10 $\mu\text{L}$  aliquot of supernatant is mixed gently by inversion with 490 $\mu\text{L}$  of toluidine blue solution in a small plastic tube. A haemocytometer chamber is prepared, the dyed cells are introduced and the liquid volume fixed by placement of a cover slip. Cells are identified into respective types and counted.

15

Cells are then diluted with physiological saline in order to obtain an average of 40,000 cells. This number may vary between ranges of 15,000-80,000 which are suitable to provide adequate responses in subsequent assays.

20

## **2. Preparing food extracts and PMA dilutions**

Freeze dried raw, whole foodstuffs are available commercially and are stored at  $-20^\circ\text{C}$  before use. To be used in the test, they are submitted to water extraction for several hours, centrifuged and filtered. The resulting filtrate constitutes the stock solution and can be used for up to one month when stored at  $2-8^\circ\text{C}$ . Appropriate aliquots of stock are then diluted with 0.9 % w/v saline to produce working solutions ready for use in the diagnostic test. These solutions are stable at  $2-8^\circ\text{C}$  for two days.

30

Other extraction procedures can be used depending to some extent on the protocols followed by the person or laboratory carrying out the test.

5 PMA (positive control) is prepared by adding 150  $\mu$ l dimethylsulfoxide (DMSO, obtained from Sigma® to 10  $\mu$ l PMA (same source) frozen stock solution (1.6 mM) in order to obtain a  $10^{-4}$  M solution. Serial 10-fold dilutions of 100  $\mu$ l PMA from  $10^{-4}$  to  $10^{-10}$  M are prepared in DMSO. Agitate on  
10 vortex for no more than 5 sec between dilutions.

### 3. Reacting the cells with PMA and a library of food extracts

15 100 $\mu$ L aliquots of a library of 81 different food extracts are dispensed into a standard 96-well microtitre plate, (polystyrene) from working dilutions prepared fresh each day from stock solutions previously prepared according to conventional methods. Twelve blanks  
20 (consisting in the suspending solution e.g. saline or saline/DMSO) are also dispensed in the plate as controls. Each of the 3 remaining wells out of the 96 available in a plate receive 1 dilution of PMA in descending order, to be chosen from  $10^{-10}$  to  $10^{-6}$  M.

25

Aliquots (100 $\mu$ L) of cell supernatant are then added to each well followed by gentle mixing for 1 min. Plates are then centrifuged at 700 rpm, for 4 min, then incubated for 1 hour at 37°C.

30

Each plate is then inverted over a sink, tapped vigorously and blotted onto tissue paper, then washed twice with saline in order to remove unbound cells.

Aliquots of 100 $\mu$ L of 1% Triton X100 detergent are then added to each well and mixed for 2 min.

5    **4. Assay for adherent PMNs.**

Since the level of LDH in TritonX-treated cell suspension is proportional to the number of cells adherent to plastic wall, the number of adherent PMNs is measured  
10 by assaying for LDH using conventional methods known per se. For example, LDH developing reagent (100 $\mu$ L) is added to each well followed by 30 min incubation, with the plates kept in darkness. The developed colour is then read in a plate reader and absorbances are recorded for a  
15 bichromatic read at 492 and 620nm. A suitable LDH reagent is that obtainable from Roche as the 'Cytotoxicity Detection Kit' Cat No 1 644 793.

Results are obtained by the following procedure:  
20 absorbance differences are calculated for each food extract between spectrophotometric reads at 492nm subtracted from the reference read at 620nm. This is termed the response. The same is repeated for the control population, (n=12), to determine a baseline response from  
25 which a significance level at 2.5 standard deviation (SD) is calculated. Reactive foods or other substances are identified as those that give a statistically significant response above the baseline response plus 2.5 SD.

30        Outlier rejection in the controls is set at 100mAU (absorbance units) and if more than 3 outliers occur then the run is rejected. A list of food extracts or other substances is compiled in ranked descending order of the



absorbance response. Foods observed to be above the level of significance are identified as those causing food intolerance. These are listed in descending order of magnitude and constitute the result. Non-reactive  
5 substances are listed separately.

### **Interpretation**

The basal state of neutrophil reactivity is indicated  
10 by the results obtained on the 12 control wells which reflect the *in vivo* state of the cells. Any ongoing inflammation (whatever the cause, reaction to foods, acute or chronic disease) or infection will result in an elevated basal reactivity.

15

The reaction to PMA indicates the cells are reactive, thus constituting an index of the quality of the blood sample and subsequent manipulations. Its level must be compared to the mean reactivity of normal neutrophils in  
20 order to detect overreactivity of the cell population upon PMA stimulation.

Reactivity to foodstuffs are analysed as described above.

25